

# **EG V1(BL)**

Self-Monitoring Blood Glucose System

## **User's Manual**



**Please read this User's Manual thoroughly  
before using your blood glucose meter.**

## **Dear EG V1(BL) SMBG System Owner,**

Thank you for using the **EG V1(BL) Self-Monitoring Blood Glucose (SMBG) System**. We designed this system to be dependable, easy-to-use, compact, lightweight and portable to help you monitor your blood glucose on a regular basis.

Please read this manual thoroughly before you begin testing. This manual provides you and your diabetes care team with important information and step-by-step direction to use the EG V1(BL) Self-Monitoring Blood Glucose System. To start testing quickly, you can also refer to the Quick Reference Guide.

Thanks again for choosing the EG V1(BL) SMBG.

## **Intended Use**

The **EG V1(BL) Self Monitoring Blood Glucose Test System** is intended for the quantitative measurement of glucose in venous whole blood or fresh capillary whole blood from fingertip, palm and forearm. Testing is done outside the body (In Vitro diagnostic use). It is indicated for use at home (over the counter [OTC]) by persons with diabetes, or in clinical settings by healthcare professionals, as an aid to monitor the effectiveness of diabetes control.

## Standard Accessories

Your new **EG V1(BL)** Blood Glucose meter and accessories work together to measure the amount of glucose in your blood. The system includes:

- **Wallet**
- **User's Manual**
- **Self-Test Log Book**
- **EG V1(BL) Blood Glucose Meter**
- **Test Strip Instructions**
- **EG Blood Glucose Test Strips (10 pcs)**
- **Code Card**
- **Quick Reference Guide**
- **CR2032 Lithium Cell 3V Battery (2 pcs)**
- **Lancets (10 pcs)**
- **Lancing device**

## Optional Accessories

- **EG Normal Control Solution**
- **EG High Control Solution**



**EG Normal Control Solution and EG High Control Solutions are available. For purchase, please contact your local dealer.**

## Why is it so important to test blood glucose regularly?

Testing your blood glucose regularly can make a big difference in how you manage your diabetes every day. We've made this SMBG system as simple as possible to help you use it regularly. Your meter is easy to use, and you can adjust the lancing device for your comfort.

## Do you need Help?

If you have questions or need assistance, please contact your healthcare professional. You can also visit [www.epsbio.com](http://www.epsbio.com) for diabetes management tools and product demonstrations.

Please complete your warranty card and mail it, so you receive the best customer service possible and news about product update.



**Although the EG V1(BL) SMBG System is easy to use, you may need to consult with your healthcare professional (this may be your doctor, pharmacist or diabetes nurse educator) for instructions on how to use the system. Only the correct use of the system will ensure accurate results.**

## **Important Information About Your New Meter**

- EG V1(BL) blood glucose meter is designed and approved for testing fresh capillary whole blood samples from your fingertip, palm and forearm. The meter is for in vitro diagnostic use ONLY (for testing outside the body). It should not be used to diagnose diabetes.
- EG V1(BL) blood glucose meter can only be used with EG Blood Glucose Test Strips. Other test strips will give inaccurate results.
- Testing is not valid for neonatal blood specimens.
- Do not disassemble the meter as this may cause damage to the components resulting in incorrect readings. Disassembling the meter will also void the warranty.
- Always keep the meter clean and store it in a safe place. Protect the meter from direct sunlight to ensure a longer lifespan.
- You should not store the meter and test strips in a car, bathroom, or refrigerator.
- Keep the meter, test strips and lancing device away from children and pets.
- You should not test critically ill patients with home-use blood glucose meters.
- Incorrect results may occur when performing the test. If you believe you are not feeling well, please contact your healthcare professional.
- Remove batteries if the meter will not be used for one month or more.

- Warning for potential biohazard: Healthcare professionals using this system on multiple patients should be aware that all products or objects that come in contact with human blood, even after cleaning, should be handled as if capable of transmitting a viral disease.
- Consult with your healthcare professional before testing on your palm or forearm.



- Do not touch the strips with wet hands.
- Do not use expired strips (the expiration date is shown on the bottle.)
- Do not bend, cut or twist the strips.
- Altitude up to 3048 meters above sea level has no effect on readings.

## Health-Related Information
















- If you are experiencing dehydration, frequent urination, low blood pressure, shock or hyperosmolar hyperglycemic nonketotic coma (HHNKC), you may get a test result that is lower than what your blood glucose really is. If you think you are dehydrated, call your doctor right away.
- If you have followed the steps in the user's manual, but still have symptoms that don't seem to match your test results, or if you have questions, please talk to your doctor or health care professional.
- Please read your test strip instruction carefully for additional health-related information.



### **Warning for potential biohazard:**









**Healthcare professionals using this system on multiple patients should handle all products or objects in contact with human blood carefully to avoid transmitting viral disease, even after cleaning.**

## Explanation of Symbols

	Consult instructions for use		Caution
	Batch code		Do not reuse
	In vitro diagnostic medical device		3V (CR2032) x 1 battery only
	Use by		Temperature limitation
	Manufacturer		Catalogue number
	Serial number		Control
	Sufficient for		Authorized representative in the European Community
	This product meets the requirements of Directive 98/79/EC in vitro diagnostic medical devices		



## Explanation of Symbols (continued)

	Discard 3 months (90 days) after opening		Large LCD screen
	Hypo alarm		Keep away from sunlight
	5-seconds result		0.6 $\mu$ L blood volume
	Alternative site testing		Human-factor flow

## Table of Content

<b>Chapter 1: Understanding Your Meter</b> .....	12
The EG V1(BL) Blood Glucose Meter .....	12
Setting The Time And Date .....	14
Setting Hypo And Hyper .....	15
Using EG Blood Glucose Test Strips .....	16
<b>Chapter 2: Control Solution Test</b> .....	17
Why Run A Control Solution Test .....	17
About The Control Solutions .....	18
Running A Control Solution Test .....	19
Understanding Control Test Results .....	21
<b>Chapter 3: Testing Your Blood Glucose</b> .....	23
Using The Lancing Device .....	23
Inserting A Lancet into the Lancing Device .....	24
Running A Blood Glucose Test With Blood From Your Fingertip .....	26
Alternative Site Testing .....	28
Running A Blood Glucose Test With Blood From Your Forearm .....	30
Running A Blood Glucose Test with Blood From Your Palm.....	31
Discarding Used Lancets .....	32
Expected Values .....	33
Unusual Test Results .....	33

Symptoms Of High Or Low Blood Glucose .....	34
Comparing Your Meter Result To Lab Results .....	35
<b>Chapter 4: Meter Memory, Setup .....</b>	<b>36</b>
Memory, Storing Test Results .....	36
Viewing Test Results .....	37
Setting The Alarm .....	38
<b>Chapter 5: Maintenance And Troubleshooting .....</b>	<b>39</b>
Inserting Batteries .....	39
Cleaning Your Meter .....	40
Maintenance And Testing .....	41
Cleaning Your Lancing Device .....	41
Screen Messages And Troubleshooting .....	42
<b>Chapter 6: Technical Information .....</b>	<b>45</b>
Specifications .....	45
Limitations.....	46
Device Information .....	48
Warranty .....	48

# Chapter 1: Understanding Your Meter

## The EG V1(BL) Blood Glucose Meter




	The Overview of the Display	Meaning
<b>A</b>	month and date	Show the moment of month and date.
<b>B</b>	Time	Show the moment of time.
<b>C</b>	Test result	Show the test result.
<b>D</b>	Error	The meter system is in fault.
<b>E</b>	Application of blood or Control Solution	Test your level of blood glucose.
<b>F</b>	Control solution test	Control Solution test, whose value won't be included in the testing memory.
<b>G</b>	Hypo	Set your hypo value. When test result is under the setting value, meter will alarm you.
<b>H</b>	Hype	Set your hype value. When test result is over the setting value, meter will alarm you.
<b>I</b>	Test records	The meter will memory automatically all of test results.
<b>J</b>	Insert test strip	Insert the test strips to start testing blood glucose.
<b>K</b>	Glucose unit	The testing unit of blood glucose is mmol/L or mg/dL.
<b>L</b>	Thermometer	The thermometer icon appears, when the temperature is too high or too low.
<b>M</b>	Batteries status	When the battery symbol appears, prepare other batteries for installation.



### The Overview of the Meter:

- A. Anti-Slippery Ring
- B. Battery cover
- C. Display
- D. Manufacturer information
- E. Serial number
- F. Strip slot
- G. Left button (◀)
- H. Right button (▶)




### The Overview of the Test Strip:


- A. Blood collection area
- B. Reaction area
- C. Holding area
- D. Strip insert direction
- E. Electrode

## Setting The Time and Date—First Time Use


Setting the current time and date in your meter is important if you use the meter memory.



1. Press “▶” until turning the meter on.  
Press “◀” once again and release to enter the setting mode.




2. The last 2-digits of the year flash at the center of the display. Press “▶” to adjust the year and Press “◀” to confirm the setting.




3. Repeat step 2 to set the date and time.  
The flashing field is the one you are currently setting.


## Setting Hypo and Hype



4. Press “▶” to adjust value of hypo and Press “◀” to confirm the setting.



5. Press “▶” to adjust value of hype and Press “◀” to confirm the setting.



6. When setting is finished, the display will be back to the testing mode of flashing strip.



**Hypo setting are from 60 to 80 mg/dL (3.3 to 4.4 mmol/L); Hype setting from 150 to 630 mg/dL (8.3 to 35.0 mmol/L).**

## Using EG Blood Glucose Test Strip

- Use only with **EG V1(BL)** Blood Glucose Meter.
- Keep the test strips in their original bottle.
- Close the bottle tightly right away, after you take a test strip out. This keeps the test strips dry.
- Use the test strip within three minutes after taking it out of the bottle.
- The test strip is for single use only. Do not reuse it.
- Record the date you open the test strip bottle. Be sure to check the expiration date on the test strip bottle. The test strip is good for three months from the date the bottle is opened or until the expiration date on the bottle, whichever comes first.
- Store the test strip bottle and your meter in a cool dry place.
- Store the test strips between 2°C ~30°C (36°F ~86°F). Do not freeze.
- Do not apply blood or control solutions to the test strip before you insert it into the meter.
- Do not touch the test strip with wet hands. Do not bend, cut, or twist the test strips.



## Chapter 2: Control Solution Testing

### Why Run A Control Solution Test

We recommend that you run the EG Normal and High Control Solutions. Such test will let you know that your meter and test strips are working properly to give reliable results. You should run the control solution tests when:

- You use the EG V1(BL) Blood Glucose Meter for the first time.
- You open a new bottle of test strips.
- You think the meter or test strips may be working incorrectly.
- You drop the meter.
- You have repeated a test and the test results are still lower or higher than expected.
- You are practicing the test procedure.
- Your test strips have been exposed to temperature outside the storage conditions (2°C to 30°C, 35.6°F to 86°F).





**Professional users are instructed to follow federal, state, and local guidelines.**

## About The Control Solutions

- Use only with EG test strips.
- Write the date you opened the bottle on the bottle label. The control solutions are good for three months from the date the bottle is opened or until the expiration date on the bottle, whichever comes first.
- Do not use control solutions that are past the expiration date?
- The control solutions can stain clothing. If you spill them, wash your clothes with soap and water.
- Close the bottle tightly after every use.
- Left over control solutions should not be added back into the control bottle.
- Store the bottles of control solutions at room temperature, between 2°C~30°C (36°F ~86°F). Do not freeze.
- If you would like to purchase EG Normal and High Control Solutions, please contact your local dealer.


## Running A Control Solution Test



1. Put a test strip into the meter in the direction of the arrow.
2. Press “▶” for 1 seconds to set “▶”.
3. Place the meter on a flat surface, like a table.
4. Remove the control solution bottle cap and wipe the tip of the bottle with a tissue.
5. Squeeze the bottle until a tiny drop forms at the tip of the bottle. Touch the drop to the Blood collection area at the end of the test strip.  
**Do not put control solution on top of the test strip.**  
When you hear the beeper, you have enough control solution in the test strip. The meter starts to count down from 5 seconds and will show the results.



6. A result appears on the display. Don't remove the test strip yet. Check if the reading falls within the range printed on the test strip bottle.



7. Remove the test strip and throw it away after you have compared the reading to the range printed on the test strip bottle.


## Understanding Control Test Results

The label on your test strip bottle shows the acceptable ranges for the Control Solutions. The result you get should be inside this range. Make sure you compare the result to the correct level of control.

When the control result is inside the range on the test strip bottle, your test strips and your meter is working properly.

If your control result is not inside the acceptable range (printed on your test strip bottle), please follow the troubleshooting checks to solve the problem.



**Control Solution values will be included in the memory and averages, if a flag icon of “” is NOT set.**

## ***Troubleshooting Checks***

- ✓ Was the test strip exposed to open air for a long period of time?
- ✓ Was the test strip bottle capped tightly?
- ✓ Was the meter functioning well?
- ✓ Are the control solutions expired or contaminated?
- ✓ Were test strips and control solutions stored in cool, dry places?
- ✓ Did you follow the testing steps properly?

## ***Action***

Repeat the control test with properly stored strips.

This will humidify strips inside. Replace the test strips.

You can use the check strip to verify the meter's functions.

Replace with a new control solution to check the performance of SMBG system.

Repeat the control test with properly stored strips or control solutions.

Read Chapter 2 "Control Solution Testing" and test again. Stop using the meter if you continue to obtain the inaccurate results. Please contact your local dealer.

## Chapter 3: Testing Your Blood Glucose

### Using The Lancing Device


- The best depth setting is the lowest number that draws enough blood for a test. Try different settings to find the one that's right for you.
- Please do not share your lancing device with anyone. And always use a new, sterile lancet. Lancets are for one time use only.





**Used test strips and lancets are considered bio-hazardous waste in accordance with local regulations and should be handled as if capable of transmitting infection. The users may discuss methods for disposing used test strips and lancets with their doctor.**

## Inserting A Lancet Into The Lancing Device


You must first load the lancet into the lancing device to get it ready for use.



1. Unscrew the Cap.



2. Insert the lancet into the lancing device firmly then twist off the protective cover.




3. Recap the front cap.




**Lancets are for use only and a new, sterile lancet should be used each time you perform a test.**






4. Select the desired penetration depth.



5. Pull on the sliding barrel of the lancing device until it clicks and then release.  
Now the lancing device is ready.  
Do not prick your finger until your meter and strip are prepared.



6. Set the lancing device aside until later in the test.



1. **Select 1-2 for soft or thin skin, 3 for average, and 4-5 for thick or calloused skin.**
2. **Lancing device and lancets are not to be shared between users. Sharing lancing devices and lancets may transmit blood borne pathogens, such as viral hepatitis.**


## Running A Blood Glucose Test With Blood From Your Fingertip



1. Wash your hands with soap and warm water. Rinse and dry thoroughly.




2. Put a test strip into the meter in the direction of the arrow. The meter turns on.



3. Swab the blood collection area on your finger with an alcohol swab. Please wait at least 5 seconds until the intended area is dry and clean before sampling the blood.
4. Place the lancing device against the pad of your finger. Press the trigger button to activate the lancing device.



**If a flag icon of “” is set, press “” for 1 seconds to cancel the setting of control solution test.**



5. When the blood drop flashes on the display, obtain a drop of blood from your fingertip.

6. Gently squeeze your finger to assist the flow of blood. This helps you get a blood drop.  
Touch the drop to the tip of the transparent window of the test strip.  
**Do not put blood on top of the strip.** Be sure to get enough blood on the strip's reaction zone.  
Otherwise, an inaccurate reading may result.

7. The meter will beep when enough blood has entered the strip's reaction zone.  
The result will appear on the display after 5 seconds.

## Alternative Site Testing (AST)

### Understanding Alternative Site Testing

#### What is AST?

Besides the fingertip, you can test the palm or forearm.

#### What is the advantage of AST?

You have the option of testing other places on your body besides the fingertip.

Consult your health care professional before you begin using the palm or forearm for testing. Blood glucose test results obtained from your palm or forearm may differ significantly from fingertip samples.

#### We strongly recommend that you:

##### Do AST ONLY in the following intervals:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

##### Do NOT use AST if:

- you think your blood glucose is low.
- you are unaware of hypoglycemia.
- your AST results do not match the way you feel.
- you are testing for hyperglycemia.
- your routine glucose results are often fluctuating.

### Fingertip test only:

- If sick
- If blood glucose is low
- After exercising
- Two hours or less after eating
- When basal insulin is most active
- After injecting rapid-acting insulin (two hours or less)
- If the blood glucose test result does not match how you feel, do a fingertip test to confirm the result again.
- Do NOT change your treatment just because of an alternate site result, do a fingertip test to confirm the result.
- If you often do not notice when your blood glucose is low, do a fingertip test.




### Caution:


- Talk with your healthcare professional before you test forearm.
- Do NOT ignore symptoms of high or low blood glucose.
- Fingertip samples are able to show the rapid change of glucose faster than forearm samples.

## Running A Blood Glucose Test With Blood From Your Forearm


Please use the clear cap with the lancing device for AST testing.



1. Massage the puncture area of forearm or forearm for a few seconds.




2. Press and hold the device with clear adjustable tip against the forearm.




3. Press the trigger button to activate the lancing device.

4. Hold the device against forearm and increase pressure until the blood sample size is sufficient.


## Running a Blood Glucose Test with Blood from Your Palm



1. Massage the puncture area of palm for a few seconds.



2. Press and hold the device with one clear adjustable tip against the palm.




3. Press the trigger button to activate the lancing device.

4. Hold the device against palm and increase pressure until the blood




**Check with your healthcare professional before testing other sites than fingertip.**


## Discarding Used Lancets



1. Unscrew and remove the cap.



2. Without touching the used lancet, stick the lancet tip into its protective cover.



3. Pointing the lancing device toward a container for sharp or biohazard material, slide the ejection button down to release the covered lancet into the container.



## Expected Values

The **EG** Blood Glucose test strips are whole-blood referenced and calibrated for easier comparison to lab results. The normal fasting blood glucose range for an adult without diabetes is 70~110 mg/dL (3.89~6.11mmol/L) \*. Two hours after meals, the blood glucose range for an adult without diabetes is less than 120 mg/dL (6.67mmol/L). For people with diabetes: please consult your doctor for the blood glucose range appropriate for you.

**Reference:** \* American Diabetes Associations: Standards of Medical Care for Patients with Diabetes Mellitus, Diabetes Care, 25(2002), p.S37

## Unusual Test Results

If your test result doesn't match the way you feel, please follow these steps:

1. Run a control test, Chapter 2, "Control Solution Testing."
2. Repeat a blood glucose test, Chapter 3, "Testing Your Blood glucose."
3. If your test results still don't reflect the way you feel, call your healthcare professional immediately.



1. **Extremely high humidity may affect the test results. A relative humidity greater than 90% may cause inaccurate results.**
2. **Hematocrit below 20% may cause higher results. Hematocrit above 60% may cause lower results.**
3. **Some studies have shown that electromagnetic fields may affect results. Do not test near an operating microwave oven.**

## Symptoms Of High Or Low Blood Glucose

Being aware of the symptoms of high or low blood glucose can help you understand your test results and decide what to do if they seem unusual. Here are the most common symptoms:

### **Greater than 13.33mmol/L(240 mg/dL) <sup>1</sup>**

#### **What It Means:**

The test result is higher than reference range.  
(3.89~6.11mmol/L or 70~110mg/dL)

#### **Symptoms:**

fatigue, increased appetite or thirst, frequent urination, blurred vision, headache, general aching, or vomiting...etc.

#### **What to Do:**

- If you are experiencing any of these symptoms, test your blood glucose.
- If the result displayed is greater than 240 mg/dL and you have symptoms of high blood glucose, contact your healthcare professional instantly.
- If the result does not match how you feel, follow the steps under "Unusual Test Results."

### **Below 3.33mmol/L(60 mg/dL) <sup>2</sup>**

#### **What It Means:**

The test result is lower than reference range.  
(3.89~6.11mmol/L or 70~110mg/dL)

#### **Symptoms:**

sweating, trembling, blurred vision, rapid heartbeat, tingling, or numbness around mouth or fingertips....etc.

#### **What to Do:**

- If you are experiencing any of these symptoms, test your blood glucose.
- If the result displayed is below 60 mg/dL and you have symptoms of low blood glucose, contact your healthcare professional instantly.
- If the result does not match how you feel, follow the steps under "Unusual Test Results."

## Comparing Your Meter Result To A Lab Result

A common question is how the blood glucose results on your meter compare to the lab results. Your blood glucose can change quickly, especially after eating, taking medication, or exercising. If you test yourself in the morning, then go to the doctor's office for a blood glucose test. The results will probably not match, even if you are fasting. This is typically not a problem with your meter, it just means that time has elapsed and your blood glucose has changed.

If you want to compare your meter result to the lab result, you must be fasting. Bring your meter to the doctor's office, and test yourself by fingertip within five minutes of having blood drawn from your arm by a healthcare professional. Keep in mind that the lab could use different technology than **EG V1(BL)** blood glucose meter, and that blood glucose meters for self testing generally read somewhat lower or higher than the lab result.

For accuracy and precision data and for important information on limitations, see the instructions that come with your test strips.

## Chapter 4: Meter Memory, Setup


### Memory, Storing Test Results

Your meter stores a maximum of 480 test results with the time and date of the test. You can review them at any time. When the memory is full, the oldest result is dropped as the newest is added, so it is very important to have the correct time and date set in the meter.



1. **Do not change your therapy based on one individual result in memory.**
2. **The memory is not lost when you replace the battery. You do need to check that the time and date are still correct. See the section "Setting the time and date" in Chapter 1.**
3. **Once 480 results are in memory, adding a new result causes the oldest one to be deleted.**


## Viewing Test Results




1. Press “▶” to view the last test result.
2. Press “◀” to view the previous results in order from the latest to the earliest.
3. Press “▶” once again to view averages of 7-day test result.
4. Press “◀” to review averages of 14-day, 30-day and 90-day test results.
5. Press “▶” to be back to the testing mode of flashing strip.

## Setting The Alarm


The meter can set up to 4 alarms, allowing you to control your test times according to your requirements.



1. Press “▶” for 5 seconds and enter the setting mode of alarm.
2. Press “◀” to set “On” for the first alarm.
3. Press “◀” to adjust time and press “▶” to confirm the setting (24-hour clock).



4. Press “▶” to confirm the first setting and then get into the second one.



5. Repeat step 1 to 4 to set other 3 alarm. If you don't want to turn on the alarm at all, turn every alarm "OFF" at the upper left corner of the display. Press “▶” to confirm the setting.

## Chapter 5: Maintenance and Troubleshooting

### Inserting Batteries



The meter uses two batteries of CR2032 Lithium Cell 3V. Batteries will normally last for over 2000 tests. Other types of CR2032 Cell 3V batteries are also acceptable, yet the capacity of test times may differ. Insert the batteries when you first use the meter or replace with new batteries when the "LP" message and the low battery symbol appear on the display.

← Low battery symbol



1. The meter won't delete earlier records after you replace batteries.
2. You should reset the time and date again after you replace the batteries.
3. CR2032 Cell 3V batteries are available in shops. You may take the old batteries with you for replacement.
4. Remove batteries when you will not be using the meter for one month or more.

## Cleaning Your Meter

Caring for your **EG V1(BL)** SMBG system does not require special cleaning. Please keep the meter free of dirt, dust, bloodstain, and water stains. Follow these guidelines carefully to help you get the best performance possible:

### Do:

- Make sure the meter is turned off.
- Gently wipe the meter's surface with a soft cloth slightly dampened.

### Do Not:

- Get any moisture in the test strip slot.
- Spray any cleaning solution directly onto the meter.
- Put the meter under water or liquid.
- Pour liquid into the meter.



## Maintenance and Testing



Your meter needs little or no maintenance with normal use. It automatically tests its own systems every time you turn it on and lets you know if something is wrong. (See "Screen Messages" and what to do about them.)

To make sure the display is working properly, turn off the meter. Press and hold "◀" or "▶" to see the complete display. All the indicators should be clear and look exactly like the picture to the left. If not, please contact your local dealer.

## Cleaning Your Lancing Device

- To clean the lancing device, wipe it with a soft cloth dampened with water and mild detergent. **DO NOT** place the entire device under water.
- To disinfect the cap after cleaning, place it in 70%~75% rubbing alcohol for 10 minutes at least once a week. Allow the cap to air-dry after disinfecting.

## Screen Messages And Troubleshooting

Never make treatment decisions based on an error message. If you have any concerns, please contact your local dealer.



### **Humidified / Used strips**

*Action:* Replace with a new strip.



### **Low power**

*Action:* Replace with new batteries.



### **System fault**

*Action:* Replace the batteries first. If 001 ERROR appears again, please contact your local dealer.



### **System fault**

*Action:* Replace the batteries first. If 005 appears again, please contact your local dealer.



The "Ht" and thermometer icon appears. Temperature is too high, outside the required range of 10°C ~40°C (50°F ~104°F). This alerts users that an incorrect result may occur if the test continues.

*Action:* Relocate the meter to a location with temperature between 10°C ~40°C (50°F ~104°F).



The "Lt" and thermometer icon appears. Temperature is too low, outside the required range of 10°C ~40°C (50°F ~104°F). This alerts users that an incorrect result may occur if the test continues.

*Action:* Relocate the meter to a location with temperature between 10°C ~40°C (50°F ~104°F).



**Err message.**  
**Not enough sample  
was drawn into the test  
strip for measurement.**  
*Action:* Discard the test  
strip and repeat the test.



**Test result is higher than  
600 mg/dL.**  
*Action:* Test again. If the  
result is still the same,  
please contact your  
healthcare professional.



**The test result is lower  
than 20 mg/dL.**  
*Action:* Test again. If the  
result is still the same,  
please contact your  
healthcare professional.

## Chapter 6: Technical Information

### Specifications

Brand name	EG V1(BL) Blood Glucose Meter	
Range	20~630 mg/dL (1.1~35.0 mmol/L)	
Response time	5 seconds	
Memory sets	480 test results	
Operating condition	Temp.	10°C ~40°C (50°F ~104°F)
	Relative Humidity	R.H. ≤ 90%
Blood sample	0.6 µL	
	Fresh blood from fingertip, palm or forearm	
Hematocrit (Hct)	20~60%	
Power	CR2032 Lithium Cell 3V (2 pcs)	
Battery life	Over 2000 tests	
Display dimension	5.1 x 3.3 mm (2.0 x 1.3 inches)	
Device dimension H × W × D	88.9 x 53.3 x 22.9 mm (3.5 x 2.1 x 0.9 inches)	
Weight	58 g (2.05 oz.) w/o batteries	
Principles	Electrochemical biosensor technology	

## Limitation

The test strips are used for fresh capillary whole blood samples.

1. DO NOT use neonate blood sample.
2. Extreme humidity may affect the results. A relative humidity greater than 90% may cause incorrect results.
3. The system should be used at a temperatures between 10°C and 40°C (50°F and 104°F). Outside this range, the system may get incorrect results.
4. DO NOT reuse the test strips. The test strips are for single use only.
5. Hematocrit: The hematocrits between 20% and 60% will not significantly affect the results. Hematocrit below 20% may cause higher results. Hematocrit above 60% may cause lower results. If you do not know your hematocrit level, please consult with your healthcare professional.
6. Altitude up to 3048 meters above sea level has no effect on readings.

Healthcare Professionals – Please note these additional Limitations

7. Patients undergoing oxygen therapy may have inaccurate results.
8. If the patient has the following conditions, the result may fail:
  - ◆ Severe dehydration
  - ◆ Severe hypotension (low blood pressure)
  - ◆ Shock
  - ◆ A state of hypoglycemic-hyperosmolar state (with or without ketosis)
9. Lipemic samples: Cholesterol level up to 500 mg/dL and triglycerides up to 1,000 mg/dL do not affect the results. Grossly lipemic patient samples have not been tested and are not recommended for testing with **EG V1(BL)** Glucose Meter.

10. Critically ill patients should not be tested with home-use blood glucose meters.
11. DO NOT use serum or plasma sample for EG V1 SMBG System.
12. DO NOT use anticoagulant NaF or potassium oxalate for venous sample preparation.
13. Interfering Substances depend on the concentration. The below substances up to the test concentration will not affect the test results.

Concentrations of the interferant tested	Bias		Glucose Level		
	80 mg/dL	250 mg/dL	500 mg/dL		
Uric acid, 8 mg/dL	2.71%	9.55%	-1.75%		
Xylose, 4 mg/dL	-5.12%	-1.64%	-4.44%		

## Device Information

---

**EG V1(BL) SMBG System,  
EG V1(BL) Blood Glucose Meter,  
EG Blood Glucose Test Strips,  
EG Normal Control Solution,  
EG High Control Solution.**



**Manufacturer:  
EPS BIO TECHNOLOGY CORP.**

No.8, R&D RD. III, Hsinchu Science Park, Hsinchu, Taiwan 30077

TEL : +886-3-6686868      FAX : +886-3-6686866

E-mail : [info@epsbio.com.tw](mailto:info@epsbio.com.tw)      Website : <http://www.epsbio.com.tw>



**Obelis s.a.**

34, Av. de Tervuren, bte 44, B-1040 Brussels, Belgium

TEL : (32) 2.732.59 54      FAX : (32) 2.732.60 03

E-mail: [mail@obelis.net](mailto:mail@obelis.net)



## Warranty

EPS warrants the original purchaser for a period of 5 years from the date of purchase. This means during the warranty period if your Self-Monitoring Blood Glucose System does not work for any reason (other than obvious abuse), EPS will replace it with a new system or an equivalent product free of charge.



## Lancing Device

Meets the requirements of MDD 93/42/EEC



**Manufacturer:**

**STERILANCE MEDICAL (SUZHOU) INC.**

No. 32 Xinlian RD. Pingjiang Suzhou P.R. China 215031

TEL : 0086 (512) 67217661 FAX : 0086 (512) 67217663

E-mail : [guopings@xinda-medical.com](mailto:guopings@xinda-medical.com)

## Lancet

Meets the requirements of MDD 93/42/EEC



**Manufacturer:**

**SAE HAN MEDICAL CORP.**

# 700-113 PUB GOT-DONG, IL SAN-GU, GOYANG-CITY, KYUNGGI-DO, KOREA

TEL : 82-31-923-4330 FAX : 82-31-923-4331

E-mail : [saehan@saehanmed.com](mailto:saehan@saehanmed.com)

## References

1. Stedman, Thomas Lathrop. Stedmans' Medical Dictionary, 27<sup>th</sup> Edition, 1999, pg. 2082.
2. American Diabetes Association, "Clinical Practice Recommendations 2003." Diabetes Care. Vol. 26, Supplement 1, pg. S22.